

APPLICATION NOTE

Performance qualification of a system measuring size and concentration of nanoparticle samples: VIDEODROP

Authors

Myriade

Anastasiia Dubrova, Marie Berger

Abstract

Viral vector applications in gene & cell therapy have been showing great potential for a variety of diseases. This involves the need for quality and safety control of the gene therapy product throughout production and quality control to ensure its progress to clinical trials. One of such regulations includes the validation of the analytical methods that are not yet standardized to be compatible with the good manufacturing practices (GMP). For instance, quantification and characterization of viral vectors production are important quality controls to perform.

Myriade is a French company that developed a new analytical device for nanoparticules characterization: Videodrop. The instrument is based on Interferometric Light Microscopy (ILM), enabling fast measurement of live concentration and size of biological nanoparticles (e.g., viral vectors) in a single drop (5-10 μ L) of a sample.

In this application note, are presented a set of tests that have been conducted on Videodrop performance for its validation as an analytical tool. The concentration and size measurements were assessed through linearity, accuracy and precision based on monodisperse samples of NIST size standard polystyrene nanobeads.

The validation procedure suggests that Videodrop is a suitable tool for **quick characterization of viral vectors**. It is an easy-to-use and fast alternative to the standard more complex and time-consuming methods.



Introduction

Cell & gene therapy has recently become one of the most rapidly evolving fields in bioengineering due to its outstanding prospective in therapeutical strategies for a huge variety of diseases. Viral vectors tend to be increasingly used as a powerful tool to introduce genes into cells ex vivo, for instance in CAR-T cell therapies. However, the success of this advanced therapeutical approach is extremely dependent on the design of standardized protocols and analytical techniques in accordance with Good Manufacturing Practice (GMP) (1).

The **control** of the production process and analytical methods used throughout manufacturing and production of viral vectors is crucial to allow final batch release. As described in EMA guidelines on the quality, non-clinical and clinical aspects of gene therapy medicinal products, the **number of particles**, the **particle size average** and **distribution** and **aggregation levels** should be determined (2).

By rapidly measuring the size and physical titer of nanoparticles in solution, Videodrop can be integrated in such quality control strategy for Drug substance (DS) and Drug Product (DP).

Validation of an analytical method for viral vector production must be done according to International Conference on Harmonization Q2 (ICH Q2) Guidelines. The validation procedure of the method performance, includes testing of various parameters like specificity, linearity, range, accuracy, precision, detection limit, quantification limit and robustness (3).


Videodrop, being one of the analytical methodologies for viral vector characterization during production, is therefore subject to such validation. In this application note, we proposed to study **linearity**, **accuracy** and **precision** of the Videodrop methodology regarding nanoparticles concentration and size. Those parameters are evaluated with size calibrated polystyrene nanoparticles.

Material and Methods

Sample preparation

The study was performed with standard size polystyrene nanobeads (3000 Series Nanosphere™

size standards, ThermoFisher, USA). 7 sizes of beads were tested (80, 100, 125, 150, 200, 300, 400 [nm]) in 7 concentrations each (8e7; 2e8; 5e8; 8e8; 2e9; 5e9; 1e10 [particles/mL]). Measurements were taken in triplicates for each sample.



Theoretical concentrations were estimated with information given by the nanospheres' supplier. Dilutions were performed in distilled water. Following suppliers' recommendations and to avoid nanoparticle aggregation, each sample was sonicated before the measurement.

Interferometric Light Microscopy (ILM)

Videodrop is a custom microscope that uses interference phenomenon to detect the light scattered by individual nanoparticles in solution.

Videos recorded by the Videodrop are processed to reveal the diffraction patterns created by the nanoparticles moving in the light path. Using this interferometric signal, nanoparticles are automatically detected and tracked to compute concentration and hydrodynamic diameter. Counting particles allows to measure the concentration, while tracking their Brownian motion allows to measure their hydrodynamical diameter.

The microscope magnification and camera speed allow to perform analysis of small sample volumes (down to 5 μ L) in less than one minute.

The Videodrop does not need any calibration or settings adjustment. It is a turnkey solution to quickly measure size and concentration of viral vector samples.

Measurement protocol

Samples were systematically diluted in distilled water to obtain the 7 studied theoretical concentrations based on the initial supplier's one. All measurements were performed to reach 300 tracked particles or 20 videos; saturation between 90 and 95%; 7 μ L sample droplets. Measurements were processed by QVIR 2.7.2 software, including likeness detector.

Linearity

Linearity regions are identified for both concentration and hydrodynamic diameter measurements based on the coefficient of determination R^2 evaluation. For size, linearity was

evaluated for 7 different bead sizes studied. The concentration linearity range was determined for 5 concentration values, excluding two of the studied concentrations (8e7 and 1e10 particles/mL). Then, the subsequent accuracy and precision analyses were conducted on these linearity regions.

Accuracy

Accuracy of measurements is addressed through the recovery rate to assess the proximity of the experimentally obtained values to the theoretical ones. It is expressed as a percentage with 100% being the ideal recovery. The recovery rate is therefore calculated as follows:

$$\text{Recovery rate} = \frac{\text{experimental}}{\text{theoretical}} \times 100\%$$

Precision



Precision study is performed to estimate the degree of scatter between a series of measurements acquired from repetitive sampling of the individual sample under the same operating conditions. Each sample, therefore, is measured 3 times. Here, the precision is expressed through a coefficient of variation (CV):

$$\text{CV} = \frac{\text{standard deviation}}{\text{mean}} \times 100\%$$

Results and Discussion

Linearity




For linearity study performance, the experimentally obtained values for both concentration and size were plotted versus the theoretical ones and the linear fit was performed.



The linear fit for concentration measurement based on 5 studied concentrations resulted in $R^2 = 0.9758$. For size measurements, the linear fit was performed on 7 studied sizes with resulting coefficient of variation $R^2 = 0.9984$. Both analyses show an excellent linearity on the range of concentration and size, both with a slope superior to 0.9.

Accuracy

Accuracy was evaluated for both size and concentration measurements plotting the recovery rate as a function of experimental results. Thus, for size measurements the ideal value or 100% recovery is represented by the supplier's indicated values. For concentration, these values correspond to theoretically derived values based on the information provided by the supplier and the subsequent successful serial dilutions. For hydrodynamic diameter accuracy study, the majority of the measurements were within 10% error, except for the two lowest sizes (80nm & 100nm) for which the error was within 20%. For concentration results, measurements are included in between $\pm 20\%$ error.




CV of concentration measurement are mostly under 20%, whereas CV of size measurement are always under 5%. Precision results are very good and promising for method validation.

Conclusion

EMA guidelines for quality control of gene therapy recommend controlling the size and number of particles of viral vectors DS, DP and critical intermediates during the process. Videodrop, measuring size and concentration of nanoparticles easily and rapidly, can naturally be integrated in **quality control strategy**. For this specific objective, analytical methods should be validated following ICH Q2 guidelines. In this application note, we proposed to determine **linearity**, **accuracy** and **precision** by analyzing polystyrene standard beads. Results of this performance study are very good and promising for future **method validation** on viral vectors for gene therapy medicinal products.

Precision

Precision was estimated through coefficient of variation (CV) based on 3 measurements of each individual sample. To represent the coefficient of variation for size and concentration measurements, the studied beads were grouped according to their size and the CV was plotted for concentration and size measurements accordingly.



	Concentration	Size
Linearity	$R^2=0.9758$	$R^2=0.9984$
Accuracy	$\pm 20\%$ error	Maj $\pm 10\%$ error
Precision	CV mostly $< 20\%$	CV $< 5\%$

References

1. The Rules Governing Medicinal Products in the European Union. In: GMP/ISO Quality Audit Manual for Healthcare Manufacturers and Their Suppliers, (Volume 2 - Regulations, Standards, and Guidelines) [Internet]. 0 éd. CRC Press; 2004 [cited 3 nov 2021]. p. 257-316. Available on: <https://www.taylorfrancis.com/books/9780203026656/chapters/10.3109/9780203026656-14>
2. Guideline on the quality, non-clinical and clinical aspects of gene therapy medicinal products. :46.
3. Q 2 (R1) Validation of Analytical Procedures: Text and Methodology. 2006;15.